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Reagent A	1 g 4,4'-dicarboxy-2,2'-biquinoline, disodium salt (Na_2 BCA; Pierce or Sigma; 26 mM final); 2 g $\text{Na}_2\text{CO}_3 \cdot \text{H}_2\text{O}$ (sodium carbonate 0.16 M final); 160 mg disodium sodium tartrate (7 mM final); 0.4 g NaOH (0.1 M final); 0.95 g NaHCO_3 (0.11 M final) brought to 100 ml with distilled water. Adjust the pH to 11.3 +/- 0.2 with 10 M NaOH. (於室溫下安定) p.s. <i>Only the disodium salt of Na_2BCA is soluble at neutral pH; the free acid is not readily soluble, even in alkali.</i>
Reagent B	0.4 g $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ (16 mM final) in 10 ml distilled water. (於室溫)
Standard working solution (SWR)	Mix 100 volumes reagent A with 2 volumes reagent B. (應為light apple-green color，於室溫下可保存一星期)

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Reagent A	8 gm sodium carbonate monohydrate, 1.6 gm sodium tartrate, with distilled water. Adjust the pH to 11.25 with 10 M NaOH
Reagent B	4 gm BCA in 100 ml distilled water.
Reagent C	0.4 gm cupric sulfate (5 x hydrated) in 10 ml water.
Working solution	Mix 1 volume reagent C with 25 volumes reagent B, then add 26 volumes reagent A to the C/B mixture. 1x C + 25x B + 26x A (順序不宜錯誤)

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The BCA assay for total protein is somewhat variable: it has differing sensitivities in response to incubation time, incubation temperature, standard protein used for calibration, and other factors (Smith et al., 1985). Certain classes of compounds such as reducing sugars and ammonium ions interfere with the assay, sometimes severely. However, if interfering compounds are eliminated (e.g., by dialysis; see Support Protocol 2), the BCA assay has a good combination of sensitivity and simplicity, and it has some advantages over the Lowry technique (see Background Information).

Reference

Stoscheck, CM. 1990. Quantitation of Protein. *Methods in Enzymology* 182: 50-69.

Smith, P.K., Krohn, R.I., Hermanson, G.T., Mallia, A.K., Gartner, F.H., Provenzano, M.D., Fujimoto, E.K., Goeke, N.M., Olson, B.J., and Klenk, D.C. 1985. Measurement of protein using bicinchoninic acid. *Anal. Biochem.* 150:76-85.

STEP 1.


Figure 4: Proteins react with alkaline copper II to produce copper I. BCA then reacts with copper I to form an intense purple color at 562 nm.

STEP 2.



Figure 3.4.3 shows a calibration plot prepared using BSA as the protein standard and a 2.1-ml final assay volume (f.a.v.). The slope of the plot is a measure of the sensitivity of the assay. The slopes of plots such as Figure 3.4.3 are also related to the absorption coefficients of the colored complexes. The units attached to the slope must take into account the extent of dilution and be consistent with the protocol used to develop the colored complexes. The value of the slope can be used to compare the different methods— e.g., Smith et al. (1985) compared the BCA and Lowry methods using BSA as calibrating protein. For the BCA method, the slope was 1.5×10^{-2} to $2.6 \times 10^{-2} A_{562} (\text{ug protein/ml f.a.v.})^{-1} \text{ cm}^{-1}$; for the Lowry method, the slope was 1.6 to $2.3 \times 10^{-2} A_{750} (\text{ug protein/ml f.a.v.})^{-1} \text{ cm}^{-1}$.

INSTRUCTIONS

Pierce BCA Protein Assay Kit

23225 23227 1296.11

Number	Description
23225	Pierce BCA Protein Assay Kit, sufficient reagents for 500 test-tube or 5000 microplate assays
23227	Pierce BCA Protein Assay Kit, sufficient reagents for 250 test-tube or 2500 microplate assays

Kit Contents:

BCA Reagent A, 1000mL (in Product No. 23225) or 500mL (in Product No. 23227), containing sodium carbonate, sodium bicarbonate, **bicinchoninic acid** and sodium tartrate in 0.1M sodium hydroxide

BCA Reagent B, 25mL, containing 4% cupric sulfate

Albumin Standard Ampules, **2mg/mL**, 10 X 1 mL ampules, containing bovine serum albumin (BSA) at 2mg/mL in 0.9% saline and 0.05% sodium azide

Storage: Upon receipt store at room temperature. Product shipped at ambient temperature.

Note: If either Reagent A or Reagent B precipitates upon shipping in cold weather or during long-term storage, dissolve precipitates by gently warming and stirring solution. Discard any kit reagent that shows discoloration or evidence of microbial contamination.

Table of Contents

Introduction	1
Preparation of Standards and Working Reagent (required for both assay procedures).....	2
Test Tube Procedure (Sample to WR ratio = 1:20)	3
Microplate Procedure (Sample to WR ratio = 1:8)	3
Troubleshooting.....	4
Related Thermo Scientific Products	5
Additional Information	6
References	8

Introduction

The Thermo Scientific™ Pierce™ BCA Protein Assay is a detergent-compatible formulation based on **bicinchoninic acid (BCA)** for the colorimetric detection and quantitation of total protein. This method combines the well-known reduction of Cu^{+2} to Cu^{+1} by protein in an alkaline medium (**the biuret reaction**) with the highly sensitive and selective colorimetric detection of the cuprous cation (Cu^{+1}) using a unique reagent containing bicinchoninic acid. The purple-colored reaction product of this assay is formed by the chelation of two molecules of BCA with one cuprous ion. This water-soluble complex exhibits a strong absorbance at 562nm that is nearly linear with increasing protein concentrations over a broad working range (20-2000 $\mu\text{g/mL}$). The BCA method is not a true end-point method; that is, the final color continues to develop. However, following incubation, the rate of continued color development is sufficiently slow to allow large numbers of samples to be assayed together.

The macromolecular structure of protein, the number of peptide bonds and the presence of four particular amino acids (cysteine, cystine, tryptophan and tyrosine) are reported to be responsible for color formation with BCA. Studies with di-, tri- and tetrapeptides suggest that the extent of color formation caused by more than the mere sum of individual color-producing functional groups. Accordingly, protein concentrations generally are determined and reported with reference to standards of a common protein such as bovine serum albumin (BSA). A series of dilutions of known concentration are prepared from the protein and assayed alongside the unknown(s) before the concentration of each unknown is determined based on the standard curve. If precise quantitation of an unknown protein is required, it is advisable to select a protein standard that is similar in quality to the unknown; for example, a bovine gamma globulin (BGG) standard (see Related Thermo Scientific Products) may be used when assaying immunoglobulin samples.

Two assay procedures are presented. Of these, the Test Tube Procedure requires a larger volume (0.1mL) of protein sample; however, because it uses a sample to working reagent ratio of 1:20 (v/v), the effect of interfering substances is minimized. The Microplate Procedure affords the sample handling ease of a microplate and requires a smaller volume (10-25 μL) of protein sample; however, because the sample to working reagent ratio is

1:8 (v/v), it offers less flexibility in overcoming interfering substance concentrations and obtaining low levels of detection.

Note: For peptide sample concentration measurements, use the Thermo Scientific™ Pierce™ Quantitative Fluorometric Peptide Assay or the Pierce™ Quantitative Colorimetric Peptide Assay (see Related Thermo Scientific Products)

Preparation of Standards and Working Reagent (required for both assay procedures)

A. Preparation of Diluted Albumin (BSA) Standards

Use Table 1 as a guide to prepare a set of protein standards. Dilute the contents of one Albumin Standard (BSA) ampule into several clean vials, preferably using the same diluent as the sample(s). Each 1mL ampule of **2mg/mL Albumin Standard** is sufficient to prepare a set of diluted standards for either working range suggested in Table 1. There will be sufficient volume for three replications of each diluted standard.

Table 1. Preparation of Diluted Albumin (BSA) Standards

Dilution Scheme for Standard Test Tube Protocol and Microplate Procedure (Working Range = 20-2,000µg/mL)			
<u>Vial</u>	<u>Volume of Diluent</u> (µL)	<u>Volume and Source of BSA</u> (µL)	<u>Final BSA Concentration</u> (µg/mL)
A	0	300 of Stock	2000
B	125	375 of Stock	1500
C	325	325 of Stock	1000
D	175	175 of vial B dilution	750
E	325	325 of vial C dilution	500
F	325	325 of vial E dilution	250
G	325	325 of vial F dilution	125
H	400	100 of vial G dilution	25
I	400	0	0 = Blank

Dilution Scheme for Enhanced Test Tube Protocol (Working Range = 5–250µg/mL)			
<u>Vial</u>	<u>Volume of Diluent</u> (µL)	<u>Volume and Source of BSA</u> (µL)	<u>Final BSA Concentration</u> (µg/mL)
A	700	100 of Stock	250
B	400	400 of vial A dilution	125
C	450	300 of vial B dilution	50
D	400	400 of vial C dilution	25
E	400	100 of vial D dilution	5
F	400	0	0 = Blank

B. Preparation of the BCA Working Reagent (WR)

- Use the following formula to determine the total volume of WR required:

$$(\# \text{ standards} + \# \text{ unknowns}) \times (\# \text{ replicates}) \times (\text{volume of WR per sample}) = \text{total volume WR required} \quad \text{Example:}$$

for the standard test-tube procedure with 3 unknowns and 2 replicates of each sample:


$$(9 \text{ standards} + 3 \text{ unknowns}) \times (2 \text{ replicates}) \times (2\text{mL}) = 48\text{mL WR required}$$

Note: 2.0mL of the WR is required for each sample in the test-tube procedure, while only 200µL of WR reagent is required for each sample in the microplate procedure.

- Prepare WR by mixing 50 parts of BCA Reagent A with 1 part of BCA Reagent B (50:1, Reagent A:B).** For the above example, combine 50mL of Reagent A with 1mL of Reagent B.

Note: When Reagent B is first added to Reagent A, turbidity is observed that quickly disappears upon mixing to yield a clear, green WR. Prepare sufficient volume of WR based on the number of samples to be assayed. The WR is stable for several days when stored in a closed container at room temperature (RT).

Procedure Summary (Test-tube Procedure, **Standard Protocol**)

Test-tube Procedure (Sample to WR ratio = 1:20)

1. Pipette 0.1 mL of each standard and unknown sample replicate into an appropriately labeled test tube.
2. Add 2.0 mL of the WR to each tube and mix well.
3. Cover and incubate tubes at selected temperature and time:
 - Standard Protocol:** 37°C for 30 minutes (working range = 20-2000 µg/mL)
 - RT Protocol:** RT for 2 hours (working range = 20-2000 µg/mL)
 - Enhanced Protocol:** 60°C for 30 minutes (working range = 5-250 µg/mL)

Notes:

- Increasing the incubation time or temperature increases the net 562nm absorbance for each test and decreases both the minimum detection level of the reagent and the working range of the protocol.
 - Use a water bath to heat tubes for either Standard (37°C incubation) or Enhanced (60°C incubation) Protocol. Using a forced-air incubator can introduce significant error in color development because of uneven heat transfer.
5. Cool all tubes to RT.
 6. With the spectrophotometer set to 562nm, zero the instrument on a cuvette filled only with water. Subsequently, measure the absorbance of all the samples within 10 minutes.
 - Note:** Because the BCA assay does not reach a true end point, color development will continue even after cooling to RT. However, because the rate of color development is low at RT, no significant error will be introduced if the 562nm absorbance measurements of all tubes are made within 10 minutes of each other.
 7. Subtract the average 562nm absorbance measurement of the Blank standard replicates from the 562nm absorbance measurement of all other individual standard and unknown sample replicates.
 8. Prepare a standard curve by plotting the average Blank-corrected 562nm measurement for each BSA standard vs. its concentration in µg/mL. Use the standard curve to determine the protein concentration of each unknown sample.

Microplate Procedure (Sample to WR ratio = 1:8)

1. Pipette 25 µL of each standard or unknown sample replicate into a microplate well (working range = 20-2000 µg/mL) (e.g., Thermo Scientific™ Pierce™ 96-Well Plates, Product No. 15041).
 - Note:** If sample size is limited, 10 µL of each unknown sample and standard can be used (sample to WR ratio = 1:20). However, the working range of the assay in this case will be limited to 125-2000 µg/mL.
2. Add 200 µL of the WR to each well and mix plate thoroughly on a plate shaker for 30 seconds.
3. Cover plate and incubate at 37°C for 30 minutes.
4. Cool plate to RT. Measure the absorbance at or near 562nm on a plate reader.

Notes:

Wavelengths from 540-590nm have been used successfully with this method. Because plate readers use a shorter light path length than cuvette spectrophotometers, the Microplate Procedure requires a greater sample to WR ratio to obtain the same sensitivity as the standard Test Tube Procedure. If higher 562nm measurements are desired, increase the incubation time to 2 hours.

Increasing the incubation time or ratio of sample volume to WR increases the net 562nm measurement for each well and lowers both the minimum detection level of the reagent and the working range of the assay. As long as all standards and unknowns are treated identically, such modifications may be useful.

5. Subtract the average 562nm absorbance measurement of the Blank standard replicates from the 562nm measurements of all other individual standard and unknown sample replicates.
6. Prepare a standard curve by plotting the average Blank-corrected 562nm measurement for each BSA standard vs. its concentration in $\mu\text{g/mL}$. Use the standard curve to determine the protein concentration of each unknown sample.

Note: If using curve-fitting algorithms associated with a microplate reader, a four-parameter (quadratic) or best-fit curve will provide more accurate results than a purely linear fit. If plotting results by hand, a point-to-point curve is preferable to a linear fit to the standard points.

Troubleshooting

Problem	Possible Cause	Solution
No color in any tubes	Sample contains a copper chelating agent	Dialyze, desalt or dilute sample Increase copper concentration in working reagent (e.g., use 50:2, Reagent A:B) Remove interfering substances from sample using Product No. 23215
Blank absorbance is OK, but standards and samples show less color than expected	Strong acid or alkaline buffer, alters working reagent pH	Dialyze, desalt, or dilute sample
	Color measured at the wrong wavelength	Measure the absorbance at 562nm
Color of samples appears darker than expected	Protein concentration is too high	Dilute sample
	Sample contains lipids or lipoproteins	Add 2% SDS to the sample to eliminate interference from lipids ³ Remove interfering substances from sample using Product No. 23215
All tubes (including blank) are dark purple	Buffer contains a reducing agent	Dialyze or dilute sample
	Buffer contains a thiol	Remove interfering substances from sample using Product No. 23215
	Buffer contains biogenic amines (catecholamines)	
Need to measure color at a different wavelength	Spectrophotometer or plate reader does not have 562nm filter	Color may be measure at any wavelength between 540nm and 590nm, although the slope of standard curve and overall assay sensitivity will be reduced

A. Interfering substances

Certain substances are known to interfere with the BCA assay including those with reducing potential, chelating agents, and strong acids or bases. Because they are known to interfere with protein estimation at even minute concentrations, avoid the following substances as components of the sample buffer:

Ascorbic acid	EGTA	Iron	Impure sucrose
Catecholamines	Impure glycerol	Lipids	Tryptophan
Creatinine	Hydrogen peroxide	Melibiose	Tyrosine
Cysteine	Hydrazides	Phenol Red	Uric acid

Other substances interfere to a lesser extent with protein estimation using the BCA assay, and these have only minor (tolerable) effects below a certain concentration in the original sample. Maximum compatible concentrations for many substances in the Standard Test Tube Protocol are listed in Table 2 (see last page of Instructions). Substances were compatible at the indicated concentration in the Standard Test Tube Protocol if the error in protein concentration estimation caused by the presence of the substance was less than or equal to 10%. The substances were tested using WR prepared immediately before each experiment. Blank-corrected 562nm absorbance measurements (for a $1000\mu\text{g/mL}$ BSA standard + substance) were compared to the net 562nm measurements of the same standard prepared in 0.9% saline. Maximum compatible concentrations will be lower in the Microplate Procedure where the sample to WR ratio is 1:8 (v/v). Furthermore, it is possible to have a substance additive affect such that even though a

single component is present at a concentration below its listed compatibility, a sample buffer containing a combination of substances could interfere with the assay.

B. Strategies for eliminating or minimizing the effects of interfering substances

The effects of interfering substances in the Pierce BCA Protein Assay may be eliminated or overcome by one of several methods.

- Remove the interfering substance by dialysis or gel filtration.
- Dilute the sample until the substance no longer interferes. This strategy is effective only if the starting protein concentration is sufficient to remain in the working range of the assay upon dilution.
- Precipitate the proteins in the sample with acetone or trichloroacetic acid (TCA). The liquid containing the substance that interfered is discarded and the protein pellet is easily solubilized in ultrapure water or directly in the alkaline BCA WR. A protocol detailing this procedure is available from our website. Alternatively, Product No. 23215 may be used (see Related Thermo Scientific Products).
- Increase the amount of copper in the WR (prepare WR as 50:2 or 50:3, Reagent A:B), which may eliminate interference by copper-chelating agents.

Note: For greatest accuracy, the protein standards must be treated identically to the sample(s).

Note: The Thermo Scientific™ Pierce™ Detergent Compatible Bradford Assay Kit (Product No. 23246) is an alternative related product compatible with a wide range of detergents.


Figure 1. Typical color response curves for BSA and BGG using the Standard Test Tube Protocol (37°C/30-minute incubation).

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
SELECTED REAGENTS FOR PRE-COLUMN DERIVATIZATION OF AMINO ACIDS

Reagent	Abbrev.	Detect.	$\lambda_{\text{abs(ex.)}}$	λ_{em}	Separation	Refs.
<i>o</i> -Phthalaldehyde	OPA	UV-FI	340/260	455	RP, MEKC	[13,52]
Phenylisothiocyanate	PTH/PTC	UV	254		RP, MEKC, CEC	[11,53-57]
4,4- <i>N,N</i> -Dimethylaminoazobenzene-4'-isothiocyanate	DABTH/DABITC	Vis	463		RP	[6]
5-Dimethylaminonaphthalene-1-sulfonyl chloride	DNS	UV-FI	360/254	570	RP, MEKC	[58-60]
4-Dimethylaminoazobenzene-4'-sulfonyl chloride	Dabsyl	Vis	420		RP	[61]
9-Fluorenylmethyl chloroformate	Fmoc	FI	265	315, 630	RP, MEKC	[62-64]
Naphthalene-2,3-dicarboxaldehyde	NDA	FI	462	490	RP, MEKC, CEC	[42,43,47,65,66]
Fluorescamine	-	UV-FI	280/390	475	RP, CE	[67,68]
Fluorescein isothiocyanate	FTH/FITC	FI	490	519	RP, MEKC	[69-71]
3-(<i>p</i> -Carboxybenzoyl)quinoline-2-carboxaldehyde	CBQ	FI	466	544	MEKC	[72]
3-(4-Carboxybenzoyl)-2-quinolinecarboxyaldehyde	CBQCA	FI	450	550	MEKC	[73]
3-(4-Tetraazolbenzoyl)-2-quinolinecarboxyaldehyde	TBQCA	Vis-FI	488/465	550	MEKC	[74]
1-Methoxycarbonylindolizine-3,5-dicarbaldehyde	IDA	Vis-FI	416/280	482	RP, CE, MEKC	[75,76]
2-(9-Anthryl)ethyl chloroformate	AEOC	UV-FI	366/256	402	RP, MEKC	[77,78]
6-Aminoquinolyl- <i>N</i> -hydroxy-succinimidyl carbamate	AQC	FI	250	400	RP, MEKC	[79,80]
Carbazole- <i>N</i> -(2-methyl)acetyl chloride	CMA-Cl	FI	335	360	RP	[81]
Acridone- <i>N</i> -acetyl chloride	ARC-Cl	FI	404	430	RP	[82]
Carbazole-9-acetyl chloride	CRA-Cl	FI	335	360	RP	[82,83]
Carbazole-9-propionyl chloride	CRP-Cl	FI	340	365	RP	[82]
Tetramethylrhodamine isothiocyanate	TRITC	FI	540	567	MEKC	[84]
2,4-Dinitrofluorobenzene	DNP	UV	254		CE	[85]
9-Isothiocyanatoacridine	-	UV	280		RP	[86]
Benzylisothiocyanate	BZITC	UV	246/238		RP	[87]

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

Mobile phase A: 1% tetrahydrofuran, 99% phosphate buffer (10 mM, pH 7.2).

Mobile phase B: methanol wavelengths of excitation 340 nm, emission 426 nm.


Time	Solvent A(%)	Solvent B(%)
0	85	15
5	80	15
15	70	30
20	70	30
25	55	45
30	50	50
45	15	85
50	10	90
51	0	100
56	0	100
57	85	15
60	85	15

Shimadzu Post-column Amino Acid Analysis System

Detection Principle

After amino acids are separated by cation exchange chromatography, aqueous sodium hypochlorite solution (to convert analytes to primary amines) and OPA/N-acetylcysteine reagent are added to the eluent, the resulting fluorescent derivatives are detected.

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Analysis Example

Two models are available - a "Na" model, capable of rapidly analyzing amino acids hydrolyzed from protein, and a "Li" model, capable of separating naturally-occurring free amino acids. The chromatograms displayed below show the separation of a standard solution of hydrolyzed protein amino acids containing 17 components, measured using the "Na" model and the "Li" model results from analyzing soy sauce diluted by factor 200, respectively.



Fig. 1. Chromatogram 2 mg/l standard solution. 1 – Aspartic acid, 2 – glutamic acid, 3 – serine, 4 – asparagine 5 – glutamine, 6 – threonine, 7 – glycine, 8 – arginine, 9 – alanine, 10 – GABA, 11 – tryosine, 12 – ethanolamine, 13 – valine, 14 – methionine, 15 – phenylalanine, 16 – histamine, 17 – isoleucine, 18 – lysine, 19 – leucine, 20 – tyramine, 21 – putrescine, 22 – cadaverine, 23 – isoamylamine, 24 – phenylethylamine, and R = reagent.

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
Analysis of Hydrolyzed Bovine Serum Albumin Using Waters AccQ•Tag™ Method for Amino Acid Analysis

Analysis of derivatized protein hydrolysate samples is simple and accurate with Waters new AccQ•Tag Method.

Conditions:

Column: Waters AccQ•Tag Column, 3.9 mm X150 mm

Column Temperature: 37°C

Mobile Phase: A ternary gradient elution profile using
Eluent A: AccQ•Tag Eluent A
Eluent B: Acetonitrile
Eluent C: MilliQ® water

Flow Rate: 1.0 ml/min

Detection: 470 Scanning
Fluorescence Detector (5µl flow cell)

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

split mode (1:100) was programmed at 3°C/min from 120 to 300°C

injector 290°C detector 300°C

W1, W2 are two tert.-butyldimethylsilyl derivatives of tryptophan

Alanine (A) Glycine (G) Valine (V) Isoleucine (I) Leucine (L) Proline (P) Methionine (M) Serine (S) Threonine (T) Phenylalanine (F) Aspartic acid (D) Glutamic acid (E) Lysine (K) Arginine (R) Histidine (H) Tyrosine (Y) Tryptophan (W) Cystine (C-C)


Equation 10: General reaction mechanism for the formation of trialkylsilyl derivatives for trimethylchlorosilane, X = Cl


Equation 19: Silylation reaction using N-methyl-N-t-butyl dimethylsilyl trifluoroacetamide (MTBSTFA) reagent: Y = O, S, NH, NR', COO, R, R' = Alk, Ar.

以上為 N,N-tert.-butyl(dimethylsilyl)trifluoroacetamide (MTBSTFA) 衍生化反應


Equation 11: Equation showing the by-product which is trimethylsilyl -acetamide from Bis(trimethylsilyl)-acetamide silylation reagent: TMS = Si(CH₃)₃, Y = O, S, NH, NR', COO, R, R' = Alk, Ar.

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Equation 12: Silylation reaction using N, N-bis(trimethyl-silyl)trifluoro-acetamide: TMS = Si(CH₃)₃, Y = O, S, NH, NR', COO, R, R' = Alk, Ar.

以上為 Bis(trimethylsilyl)trifluoroacetamide (BSTFA) 衍生化反應

Equation 15: Silylation reaction using Hexamethyldisilzane (HMDS): Y = O, S, NH, NR', COO, R, R' = Alk, Ar.

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

FIG. 2. An example chromatogram obtained from a standard solution containing a mixture of compounds at 45 μmol. (1) Lactate, (2) pyruvate, (3) alanine, (4) glycine, (5) valine, (6) leucine, (7) oxalate, (8) isoleucine, (9) malate, (10) proline, (11) methionine, (12) serine, (13) threonine, (14) phenylalanine, (15) asparagine, (16) hydroproline, (17) cysteine, (18) ornithine, (19) lysine, (20) glutamine, (21) arginine, (22) ascorbate, (23) histidine, (24) tryptophan. Other conditions as described in the legend to Fig. 1.

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Amylase (DU) “Dextrinizing unit” -- One FCC alpha-amylase Dextrinizing Unit (DU) is defined as the quantity of alpha-amylase that will dextrinize soluble starch in the presence of an excess of beta-amylase at the rate of one gram per hour at 30°C. The assay is based on the time required for an alpha-amylase to hydrolyze a limit dextrin substrate to a defined blue value. The degree of hydrolysis is determined by comparing the blue value of the hydrolyzed limit dextrin substrate with that of a color standard.

0.4mL KP-buffer (0.1 M pH 6.9) + 0.2 mL starch (2 mg/mL, in boiling water 10 min then cooling) + 0.2mL Enzyme \rightarrow 50°C (or 37°C)、10 min \rightarrow 0.2mL 1N HCl 終止反應

\rightarrow +1.0 mL iodine solution (5mM iodine+ 5mM KI) \rightarrow 測 580 nm (有沉澱時，離心後再測)


Standard curve (0-0.4 mg)

Starch (mg)	0	0.1	0.2	0.3	0.4	
starch (2 mg/mL)	0 μ L	50 μ L	100 μ L	150 μ L	200 μ L	
KP-buffer (μ L)	800 μ L	750 μ L	700 μ L	650 μ L	600 μ L	

0.8 mL Standard solution + 0.2mL 1N HCl + 1.0 mL iodine solution (5mM iodine+ 5mM KI)

→ 測 580 nm

One Dextrinizing unit (DU)= 1mg starch reduce/min-g

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*酵素活性單位定義:1U= liberate 1.0 mg of maltose from starch /min (at pH 6.9 at 50°C)

1U= liberate 1.0 μ mole of maltose from starch /min (at pH 6.9 at 50°C)

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
1U= liberate 1.0 mg of maltose from starch per min at 50°C and pH 6.9

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